




ANTECH IMAGING SERVICES

<https://info.antechimaging.com>

## How to submit a PennHIP consult (Study stored in PACS)

1. Log in to your account on <https://info.antechimaging.com> by entering your Username and Password and click “Sign In.”

2. Study Date range is defaulted to the last 7 days. You can change the “From Study Date” and “To Study Date” to your desired search dates. You can also search by any of the other boxes. Then click “Go” to show your studies.

3. Find the patient for which you want to send a consult request. Click the  icon or choose “Submit a Consult” from the Select Action drop-down list.

4. Ensure the patient information on the next page is correct and edit if needed. Fields marked with \* are required. **Note:** Patient ID is the patient ID your hospital uses and will ensure previous/future studies are linked to that patient in our system. Select “PennHIP” then click “Submit” to go to the next page.



ANTECH IMAGING SERVICES

<https://info.antechimaging.com>





## How to submit a PennHIP consult (Study stored in PACS)

5. Complete the required fields of the PennHIP form (indicated by a \*). If there is additional information to convey, this can be typed into the “Comments” box, directly under the Restraint Drugs box. Click **“Submit with selected radiographs.”**

Please Complete All Sections and Submit with Radiographs * - Required			
Consultation Type		Current Consult Status: Pending	
Submission for PennHIP Certification			
PennHIP Evaluation Information - AIS Demonstration Center			
Member Name*: Dr. Joe Demonstration - Member# CERT - Distractor#			
Radiograph Date*: 02/03/2024 (MM/DD/YYYY)	Patient Weight*: <input type="text"/>	Patient ID*: 239095932	
Clinical Signs*: <input type="radio"/> Yes <input type="radio"/> No <input type="radio"/> Not Evaluated	Restraint Drugs*: <input type="text"/>		
Severity: <input type="radio"/> Mild <input type="radio"/> Moderate <input type="radio"/> Severe (See Chapter 7 in Manual)	Comments: <input type="text"/>		
Owner Information			
Last Name*: Patricia		First Name*: Chestnut	
Patient Information			
Registered Name*: <input type="text"/>		Patient Name*: Amerie	
Breed*: <input type="text"/> Please Select A Breed	Gender*: <input type="text"/> Female	Date Of Birth*: 02/03/2021 (MM/DD/YYYY)	
Tattoo Number*: <input type="text"/>	Microchip Number*: <input type="text"/>		
Registration Number*: <input type="text"/>	Sire Registration Number*: <input type="text"/>	Dam Registration Number*: <input type="text"/>	
Registry/Kennel Club*: <input type="text"/>	* If your Registry/Kennel Club is not listed, contact 877-727-6800		
PennHIP ID: <input type="text"/>			
IMPORTANT: Has this dog had hip surgery? <input type="radio"/> Yes <input type="radio"/> No <input type="radio"/> Unknown			
Has THIS dog had a PennHIP radiograph before? <input type="radio"/> Yes <input type="radio"/> No <input type="radio"/> Unknown			If yes, When: <input type="text"/>
Has This Dog Suffered Hip Trauma? <input type="radio"/> Yes <input type="radio"/> No <input type="radio"/> Unknown			Comments: <input type="text"/>
Submit with selected radiographs			
Save Changes			
Cancel This Consultation			

8. A Consultation Confirmation will appear. If desired, you may print out this page for your records to send with the images. You are now done and may click on **“Back to PACS”** at the bottom to go back to your PACS Workbench.

**Note:** Check report status on your PACS workbench.

-  -The images are stored, but have NOT been sent for a consult, or the case is in a Returned state.
-  -The consult is submitted, but no (or not all) images have been received.
-  -The case is submitted and is “In Progress” of being read.
-  -There is a final Report. Click this icon to display your report in a separate pop-up window.