

Maintenance and Control Log

Control Lot# ____ Date Opened ____ Control Lot# ___ Date Opened ____ Date Opened ____

Month/Year:	

					control coth Date opened	
Day	Run Background	Run Control	Run Reticulocyte Control	Probe Cleanser*	Comments	Initials
1	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
2	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail			
3	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
4	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
5	□ Pass □ Fail	□ Pass □ Fail	□ Pass □ Fail			
6	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
7	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail			
8	□ Pass □ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail			
9	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail			
10	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail			
11	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail			
12	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
13	□ Pass □ Fail	□ Pass □ Fail	□ Pass □ Fail			
14	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
15	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
16	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
17	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail	_		
18	☐ Pass ☐ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
19	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐ Pass ☐ Fail			
20	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
21	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
22	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
23	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
24	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
25	□ Pass □ Fail	□ Pass □ Fail	☐ Pass ☐ Fail			
26	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail	_		
27	□ Pass □ Fail	□ Pass □ Fail	□ Pass □ Fail			
28	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
29	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
30	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			
31	☐ Pass ☐ Fail	□ Pass □ Fail	□ Pass □ Fail			

Heska recommends performing quality control daily. Store controls at refrigerated temperature. Discard if frozen. Bring controls to room temperature and mix well prior to use. Return controls to the refrigerator immediately after use. Hematology control is stable for 14 days after opening.

*Probe Cleanser is recommended every 7 days after installation.



For further assistance, please call Heska's Customer Support Services

US 800 464 3752 www.heska.com CA 866 382 6937 www.heskavet.ca AU 1300 437 522 www.heska.com.au